

TUTORIAL 8 KRAKEN

BY CODY SCULLY & CHRISTOPHER UZOKWE

TAXONOMIC CLASSIFICATION




Blast (Magablast)

The Naïve Bayes Classifier (NBC)

PhymmBL


MatePhlAn

Kraken



OVERVIEW OF KRAKEN

- This is a Program that reads in Group of unmapped data
- Classifies the data using k-mers to map the LCA
- Creates a Taxonomy tree with the Classifications
- Removes any unclassified leaves from the tree



*
<G.<G<AGGII...

- sabre
- fastx_dem
- idemp
- fastx barco

read-based analysis

no-assembly path

Some tools:

- TIPP/SEPP
- metaphlan2
- humann2
- sourmash
- kraken

Analysis

Some tools:

- phyloseq
- Breakaway
- DivNet
- CORNCOB
- SpiecEasi
- MaAsLin
- DESeq2

Count Table

	Sample_A	Sample_B	...
obj_1	0	428	...
obj_2	306	323	...
obj_3	217	1	...
...

MetaQUAST is a great tool for comparing assemblies

assembly path

(co)-as

Some assemblers:

- Megahit (assembly)
- SPAdes (assembly)
- idba-ud (assembly)
- MetAMOS (assembly)
- MetaCompas
- Metagenome


Gene calling
Functional/taxonomic
profiling

OVERVIEW OF METAGENOMICS WORKFLOW

QUERY SEQUENCE


- The First step that Kraken takes is to sequence the data set
 - This is done by a library of metagenome data
- Every k-Mer is either
 - **Classified** with the data from the library
 - **Unclassified** if nothing is found for them.

Query sequence



TAXONOMY TREE


- Each K-mer is mapped to the lowest common ancestor (LCA) of the genomes that contains that k-mer



Taxonomy tree


CLASSIFICATION TREE AND PATH

- Finally the unclassified K-mers are removed and the tree is reduced and weighted/scored



KRAKEN VARIANCE

- Kraken
 - Needs 70 GB of RAM
- Kraken-Q/MiniKraken-Q
 - Faster then Kraken /MiniKraken
 - Decreased Accuracy
- MiniKraken
 - Reduced database
 - Needs 4 GB of RAM
 - Less overall Sensitivity (~11-25%)
- Kraken-GB
 - Uses GenBank Database
 - Very high sensitivity





KRAKEN2 TUTORIAL



INSTALLATION & DOWNLOADS

- Download Kraken2 Docker:

```
$ conda create --yes -n kraken kraken2 bracken  
$ conda activate kraken
```

- (Or) find it on Picotte:

```
biobakery_workflows.sif  edirect_latest.sif  metaspades_latest.sif  nextflow  qiime2_latest.sif  
conda-qiime2_latest.sif  kraken2_latest.sif  miniconda3  nfcore-magbusco-1.2.0.img  README  
diamond_latest.sif      metabat_latest.sif  modulefiles  qiime  
[cnu25@picotte001 containers]$
```

- Must also download taxo database:

```
$ curl -O ftp://ftp.ccb.jhu.edu/pub/data/kraken2_dbs/minikraken2_v2_8GB_201904_UPDATE.tgz  
  
# alternatively we can use wget  
$ wget ftp://ftp.ccb.jhu.edu/pub/data/kraken2_dbs/minikraken2_v2_8GB_201904_UPDATE.tgz
```

- Also available on picotte

```
k2_pluspf_20210127  k2_pluspf_20210127.tar.gz  
[cnu25@picotte001 kraken2_db_plus_protozoa_fungi]$
```

RUN KRAKEN2 CONTAINER IN SINGULARITY

- Request node to run on:

```
[cnu25@picotte001 data]$ srun --nodes=1 --ntasks=1 --cpus-per-task=32 --mem=120GB --time=00:30:00 --pty /bin/bash
```

- Run the docker using singularity, bind to our groups folder

```
[cnu25@node003 eces450650Grp]$ singularity run --bind /ifs/groups/eces450650Grp/ containers/kraken2_latest.sif
```

- We are using the run command to execute through the kraken container

```
Usage:
singularity [global options...]

Description:
Singularity containers provide an application virtualization layer enabling
mobility of compute via both application and environment portability. With
Singularity one is capable of building a root file system that runs on any
other Linux system where Singularity is installed.

Options:
-c, --config string    specify a configuration file (for root or
                        unprivileged installation only) (default
                        "/etc/singularity/singularity.conf")
-d, --debug            print debugging information (highest verbosity)
-h, --help             help for singularity
                        --nocolor    print without color output (default False)
-q, --quiet            suppress normal output
-s, --silent           only print errors
-v, --verbose          print additional information

Available Commands:
build      Build a Singularity image
cache      Manage the local cache
capability Manage Linux capabilities for users and groups
config     Manage various singularity configuration (root user only)
delete     Deletes requested image from the library
exec       Run a command within a container
help       Help about any command
inspect    Show metadata for an image
instance   Manage containers running as services
key         Manage OpenPGP keys
oci         Manage OCI containers
plugin     Manage Singularity plugins
pull       Pull an image from a URI
push       Upload image to the provided URI
remote     Manage singularity remote endpoints, keyserver and OCI/Docker registry credentials
run         Run the user-defined default command within a container
run-help   Show the user-defined help for an image
search     Search a Container Library for images
shell      Run a shell within a container
sif         siftool is a program for Singularity Image Format (SIF) file manipulation
sign       Attach digital signature(s) to an image
test       Run the user-defined tests within a container
verify     Verify cryptographic signatures attached to an image
version    Show the version for Singularity

Examples:
$ singularity help <command> [<subcommand>]
$ singularity help build
$ singularity help instance start

For additional help or support, please visit https://www.sylabs.io/docs/
[cnu25@node003 eces450650Grp]$
```

USAGE NOTES (KRAKEN)

- `--use-names`: print scientific names
- `--db`: Database reference file path
- `--fastq-input` [deprecated] : We are using fastq formatted files
- `--paired`: We are dealing with paired end data
- `--report FILE`: This provides us with a sample-wide report

```
Singularity> kraken2 help
kraken2: Must specify DB with either --db or $KRAKEN2_DEFAULT_DB
Singularity> kraken2 --help
Usage: kraken2 [options] <filename(s)>

Options:
  --db NAME                Name for Kraken 2 DB
                           (default: none)
  --threads NUM            Number of threads (default: 1)
  --quick                  Quick operation (use first hit or hits)
  --unclassified-out FILENAME
                           Print unclassified sequences to filename
  --classified-out FILENAME
                           Print classified sequences to filename
  --output FILENAME        Print output to filename (default: stdout); "-" will
                           suppress normal output
  --confidence FLOAT       Confidence score threshold (default: 0.0); must be
                           in [0, 1].
  --minimum-base-quality NUM
                           Minimum base quality used in classification (def: 0,
                           only effective with FASTQ input).
  --report FILENAME        Print a report with aggregate counts/clade to file
  --use-mpa-style          With --report, format report output like Kraken 1's
                           kraken-mpa-report
  --report-zero-counts    With --report, report counts for ALL taxa, even if
                           counts are zero
  --report-minimizer-data  With --report, report minimizer and distinct minimizer
                           count information in addition to normal Kraken report
  --memory-mapping        Avoids loading database into RAM
  --paired                The filenames provided have paired-end reads
  --use-names             Print scientific names instead of just taxids
  --gzip-compressed       Input files are compressed with gzip
  --bzip2-compressed       Input files are compressed with bzip2
  --minimum-hit-groups NUM
                           Minimum number of hit groups (overlapping k-mers
                           sharing the same minimizer) needed to make a call
                           (default: 2)
  --help                  Print this message
  --version               Print version information

If none of the *-compressed flags are specified, and the filename provided
is a regular file, automatic format detection is attempted.
```

KRAKEN2 COMMAND, I/O, RUNTIME

```
Singularity> kraken2 --use-names --threads 4 --db ./data/kraken2_db_plus_protozoa_fungi/k2_pluspf_20210127 --fastq-input -report ./ECES450650_SP21/Tutorial8/reporttest --paired ./data/mappings/evol1.sorted.unmapped.R1.fastq ./data/mappings/evol1.sorted.unmapped.R2.fastq> ./ECES450650_SP21/Tutorial8/evol1.krakentest
```

- Input files: /ifs/groups/eces450650Grp/data/mappings/evol1.sorted.unmapped.R1.fastq
/ifs/groups/eces450650Grp/data/mappings/evol1.sorted.unmapped.R1.fastq
- Output files: ./ECES450650_SP21/Tutorial8/evol1.kraken ← classifications
./ECES450650_SP21/Tutorial8/report ← report

Execution time:

```
17692 sequences (0.85 Mbp) processed in 0.039s (26969.5 Kseq/m, 1291.14 Mbp/m)  
877 sequences classified (4.96%)  
16815 sequences unclassified (95.04%)
```

RESULTS & UNDERSTANDING - CLASSIFICATIONS

1.C/U: one letter code indicating that the sequence was either classified or unclassified.

2.The sequence ID, obtained from the FASTA/FASTQ header.

3.The taxonomy ID Kraken2 used to label the sequence; this is **0** if the sequence is unclassified and otherwise should be the NCBI Taxonomy identifier.

4.The length of the sequence in bp.

5.A space-delimited list indicating the lowest common ancestor (in the taxonomic tree) mapping of each k-mer in the sequence. For example, **562:13 561:4 A:31 0:1 562:3** would indicate that:

- the first 13 k-mers mapped to taxonomy ID #562
- the next 4 k-mers mapped to taxonomy ID #561
- the next 31 k-mers contained an ambiguous nucleotide
- the next k-mer was not in the database
- the last 3 k-mers mapped to taxonomy ID #562

Source:

<https://genomics.readthedocs.io/en/latest/ngs-taxonomic-investigation/index.html#wood2014>

```
U NS500207:12:H04WYAFXX:4:11509:6532:2207 unclassified (taxid 0) 16|19 |:|
U NS500207:12:H04WYAFXX:2:11108:4165:7520 unclassified (taxid 0) 21|17 |:|
U NS500207:12:H04WYAFXX:4:21512:24440:6715 unclassified (taxid 0) 17|16 |:|
U NS500207:12:H04WYAFXX:4:21511:14152:1702 unclassified (taxid 0) 15|23 |:|
U NS500207:12:H04WYAFXX:1:11107:18692:15804 unclassified (taxid 0) 16|15 |:|
U NS500207:12:H04WYAFXX:4:11508:6821:7156 unclassified (taxid 0) 17|16 |:|
U NS500207:12:H04WYAFXX:4:21501:17810:20140 unclassified (taxid 0) 17|18 |:|
U NS500207:12:H04WYAFXX:4:21508:12576:14696 unclassified (taxid 0) 17|17 |:|
U NS500207:12:H04WYAFXX:3:11602:22795:11585 unclassified (taxid 0) 15|17 |:|
U NS500207:12:H04WYAFXX:4:21509:18909:3154 unclassified (taxid 0) 16|22 |:|
U NS500207:12:H04WYAFXX:1:21109:21187:4328 unclassified (taxid 0) 15|16 |:|
U NS500207:12:H04WYAFXX:4:11409:22838:17093 unclassified (taxid 0) 18|20 |:|
U NS500207:12:H04WYAFXX:2:21103:12426:6694 unclassified (taxid 0) 24|19 |:|
U NS500207:12:H04WYAFXX:4:21510:14090:16787 unclassified (taxid 0) 16|16 |:|
U NS500207:12:H04WYAFXX:4:21508:1625:10321 unclassified (taxid 0) 16|16 |:|
U NS500207:12:H04WYAFXX:2:11308:1672:19011 unclassified (taxid 0) 16|15 |:|
U NS500207:12:H04WYAFXX:4:11504:25226:9766 unclassified (taxid 0) 15|17 |:|
U NS500207:12:H04WYAFXX:4:21407:22845:18028 unclassified (taxid 0) 17|15 |:|
U NS500207:12:H04WYAFXX:3:11510:6276:5073 unclassified (taxid 0) 17|27 |:|
U NS500207:12:H04WYAFXX:4:21509:15574:13445 unclassified (taxid 0) 17|17 |:|
U NS500207:12:H04WYAFXX:1:21310:2425:18928 unclassified (taxid 0) 20|15 |:|
U NS500207:12:H04WYAFXX:3:21402:24567:16422 unclassified (taxid 0) 15|16 |:|
U NS500207:12:H04WYAFXX:1:11211:25220:12773 unclassified (taxid 0) 15|17 |:|
U NS500207:12:H04WYAFXX:4:21503:7033:1837 unclassified (taxid 0) 15|15 |:|
U NS500207:12:H04WYAFXX:4:11608:21357:5621 unclassified (taxid 0) 19|17 |:|
U NS500207:12:H04WYAFXX:3:21405:13435:4254 unclassified (taxid 0) 18|16 |:|
U NS500207:12:H04WYAFXX:4:21512:15141:3807 unclassified (taxid 0) 16|24 |:|
U NS500207:12:H04WYAFXX:4:21503:14233:5503 unclassified (taxid 0) 17|16 |:|
U NS500207:12:H04WYAFXX:4:21501:1687:12407 unclassified (taxid 0) 16|16 |:|
U NS500207:12:H04WYAFXX:2:21311:19996:18525 unclassified (taxid 0) 18|15 |:|
U NS500207:12:H04WYAFXX:2:21306:11783:14063 unclassified (taxid 0) 18|20 |:|
U NS500207:12:H04WYAFXX:4:21512:9336:3276 unclassified (taxid 0) 17|15 |:|
U NS500207:12:H04WYAFXX:1:21105:6962:1196 unclassified (taxid 0) 24|18 |:|
U NS500207:12:H04WYAFXX:4:21508:21403:13767 unclassified (taxid 0) 16|23 |:|
U NS500207:12:H04WYAFXX:1:21310:13009:11906 unclassified (taxid 0) 15|17 |:|
U NS500207:12:H04WYAFXX:1:11108:12003:16084 unclassified (taxid 0) 18|22 |:|
U NS500207:12:H04WYAFXX:4:11504:21966:11645 unclassified (taxid 0) 17|15 |:|
U NS500207:12:H04WYAFXX:4:21511:1359:5478 unclassified (taxid 0) 17|16 |:|
U NS500207:12:H04WYAFXX:1:21301:19565:6703 unclassified (taxid 0) 27|17 |:|
U NS500207:12:H04WYAFXX:4:21502:11293:19099 unclassified (taxid 0) 17|16 |:|
U NS500207:12:H04WYAFXX:3:11506:20048:6180 unclassified (taxid 0) 16|16 |:|
U NS500207:12:H04WYAFXX:3:21608:10117:19112 unclassified (taxid 0) 24|27 |:|
U NS500207:12:H04WYAFXX:4:11505:1634:5673 unclassified (taxid 0) 17|15 |:|
U NS500207:12:H04WYAFXX:3:21502:9704:14065 unclassified (taxid 0) 18|17 |:|
U NS500207:12:H04WYAFXX:4:11508:26200:14533 unclassified (taxid 0) 17|17 |:|
```


RESULTS & UNDERSTANDING – REPORT SUMMARY

1. **Percentage** of reads covered by the clade rooted at this taxon
2. **Number of reads** covered by the clade rooted at this taxon
3. **Number of reads** assigned directly to this taxon
4. A rank code, indicating **(U)nclassified, (D)omain, (K)ingdom, (P)hylum, (C)lass, (O)rder, (F)amily, (G)enus, or (S)pecies**. All other ranks are simply “-”.
5. NCBI Taxonomy ID
6. The indented scientific name

Source:

<https://genomics.readthedocs.io/en/latest/ngs-taxonomic-investigation/index.html#wood2014>

95.04	16815	16815	U	0	unclassified
4.96	877	1	R	1	root
4.93	872	0	R1	131567	cellular organisms
4.54	804	0	D	2	Bacteria
4.37	774	0	P	1224	Proteobacteria
4.35	769	0	C	1236	Gammaproteobacteria
4.31	763	0	O	91347	Enterobacteriales
4.31	762	726	F	543	Enterobacteriaceae
0.14	25	0	G	561	Escherichia
0.11	20	19	S	562	Escherichia coli
0.01	1	1	S1	2027293	Escherichia coli O8:H8
0.02	4	0	G1	2608889	unclassified Escherichia
0.02	4	4	S	2044467	Escherichia sp. E4742
0.01	1	1	S	1499973	Escherichia marmotae
0.03	5	2	G	620	Shigella
0.01	2	1	S	623	Shigella flexneri
0.01	1	1	S1	374923	Shigella flexneri 1a
0.01	1	1	S	622	Shigella dysenteriae
0.01	2	1	G	413496	Cronobacter
0.01	1	0	G1	2649764	unclassified Cronobacter
0.01	1	1	S	1906275	Cronobacter sp. JZ38
0.01	1	0	G	579	Kluyvera
0.01	1	1	S	61648	Kluyvera intermedia
0.01	1	1	G	590	Salmonella
0.01	1	0	G	158483	Cedecea
0.01	1	1	S	158822	Cedecea neteri
0.01	1	0	G	1330547	Kosakonia
0.01	1	0	G1	2632876	unclassified Kosakonia
0.01	1	1	S	2725560	Kosakonia sp. SMBL-WEM22
0.01	1	0	F	1903410	Pectobacteriaceae
0.01	1	0	G	122277	Pectobacterium
0.01	1	1	S	1905730	Pectobacterium parmentieri
0.03	6	0	O	72274	Pseudomonadales
0.03	6	0	F	468	Moraxellaceae
0.03	6	2	G	469	Acinetobacter
0.02	4	4	S	40215	Acinetobacter junii
0.02	3	0	C	28216	Betaproteobacteria
0.02	3	0	O	80840	Burkholderiales
0.01	2	0	F	80864	Comamonadaceae
0.01	2	0	G	80865	Delftia
0.01	1	1	S	80866	Delftia acidovorans
0.01	1	0	G1	2613839	unclassified Delftia
0.01	1	1	S	1920191	Delftia sp. HK171
0.01	1	0	F	119060	Burkholderiaceae
0.01	1	0	G	48736	Ralstonia

SOURCES

- “6. Taxonomic investigation¶,” 6. *Taxonomic investigation - Genomics Tutorial 2020.2.0 documentation*. [Online]. Available: <https://genomics.readthedocs.io/en/latest/ngs-taxonomic-investigation/index.html#kraken2>. [Accessed: 20-May-2021].
- D. E. Wood and S. L. Salzberg, “Kraken: ultrafast metagenomic sequence classification using exact alignments,” *Genome Biology*, 03-Mar-2014. [Online]. Available: <https://genomebiology.biomedcentral.com/articles/10.1186/gb-2014-15-3-r46>. [Accessed: 20-May-2021].
- D. Wood, “DerrickWood/kraken2,” *GitHub*. [Online]. Available: <https://github.com/DerrickWood/kraken2>. [Accessed: 20-May-2021].
- M. Lee, “Genomics,” *Happy Belly Bioinformatics*, 2021. [Online]. Available: <https://astrobiomike.github.io/genomics/>. [Accessed: 20-May-2021].